# theory\_gaussian

May 21, 2020

## 1 Theory: Bessel-Gauss beam

[4]: # from beams package ./beams/test\_gaussian.py
from test\_gauss import plot\_pupil\_function, plot\_gaussian

#### 1.0.1 References:

wiki


#### 1.1 The Gaussian beam

The Gaussian beam is described at the focus by the Gaussian function:

$$E(x) = \exp\left(\frac{-x^2}{w_0^2}\right) \,,$$

where  $w_0$  is the beam waist radius.

Note, at  $x = w_0$ ,  $E(w_0) = \exp(-1) = 1/e$ ; therefore,  $w_0$  describes the radius from centre where the amplitude drops to 1/e of its maximum. Also note that intensity (which we measure on the camera)  $I = |E|^2$ ; thus,  $I(w_0) = 1/e^2$ , and  $w_0$  describes the  $1/e^2$  radius. Finally,  $2w_0$  is the  $1/e^2$  diameter, and  $2w_0 = FWHM \times \sqrt{2}/\ln(2)$ .



### 1.1.1 Gaussian beam in the pupil plane

Let's assume that the Gaussian beam in the pupil plane is given as:

$$E_2(x) = \exp\left(\frac{-x^2}{w_0^2}\right) .$$

## 1.1.2 Gaussian beam in the image plane

Using the Fourier transforming property of the lens, we can recover the amplitude at the image plane. Ignoring the amplitude component:

$$E_1(u) = \exp(-\pi^2 w_0^2 u^2),$$

which is evaluated at  $u = x_1/(f\lambda)$ . Thus, the field at the focus becomes

$$E_1(x_1) = \exp(-\pi^2 w_0^2 x_1^2 / (f\lambda)^2) = \exp(-x_1^2 / w_1^2),$$

where  $w_1 = f\lambda/(\pi w_0)$ .

#### 1.1.3 Gaussian beam propagation at the sample

Importantly, we know that the beam waist at focus (1/e radius of amplitude) is  $w_1 = f\lambda/(\pi w_0)$ .

We can figure out beam propagation from the Fresnel diffraction solution to the Helmholtz equation, but this is done for us and is well-known for a Gaussian beam. The Gaussian beam amplitude is given as:

$$E_1(x,z) = \frac{w_1}{w(z)} \exp\left(\frac{-x^2}{w(z)^2}\right) \exp\left(-i(kz + k\frac{x^2}{2R(z)} - \psi(z))\right),$$


where  $w(z) = w_1 \sqrt{1 + (z/z_R)^2}$ ,  $z_R = \pi w_1^2 n/\lambda$  is the Reyleigh range,  $R(z) = z(1 + (z_R/z)^2)$  is the radius of curvature, and  $\psi(z) = \arctan(z/z_R)$  is the Guoy phase.

## [5]: plot\_gaussian()

### 2.174693142407658e-06

F:\Work\Projects\deep-learning\deep-learning-lsm\beams\test\_gauss.py:60:
RuntimeWarning: divide by zero encountered in double\_scalars
rz = zv\*(1+(zr/zv)\*\*2)

F:\Work\Projects\deep-learning\deep-learning-lsm\beams\test\_gauss.py:60:
RuntimeWarning: invalid value encountered in double\_scalars
rz = zv\*(1+(zr/zv)\*\*2)



### 1.1.4 Implementation note

Note that in the formula for  $E_1$ , if z = 0, R(z) is a division by 0.

We make a modified parameter  $zR(z)=(z^2+z_R^2)$ , and modify the last exponent to include  $kx^2z/(2(zR_z(z)))$  which evaluates to 0 when z=0.